




**KHIMOD**  
ALCEN

**CONTINUOUS  
HYDROGENATION  
CASES STUDIES**



# CP2M Multiphase Catalytic Millifluidic Platform


## Reactor Characterization

- Heat transfer
- Residence time distribution
- Gas/Liquid & Liquid /Solid mass transfer


## Heterogeneous Catalysis

- Batch to flow
- Process intensification
- Reactive Liquid/Liquid extraction

### Pulsed experiments



### Transient experiments



### Axel'One Campus



### CPE-Lyon



➤ These case studies have been conducted in the frame of a collaboration between KHIMOD and the laboratory CP2M



# Definition and importance of heterogeneous hydrogenation catalysis

**Petroleum Refining:** Hydrocracking ; Desulfurization

**Chemical Manufacturing:** Ammonia Synthesis; Cyclohexane Production

**Fine Chemicals / Intermediates:** Hydrogenation of Nitro; Synthesis of Alcohols


**Food Industry:** Hydrogenation of Oils

**Pharmaceutical Industry:** Synthesis of Active Pharmaceutical Ingredients (APIs)

Key benefits of flow hydrogenation :

- Scalability (numbering-up)
- More precise control over reaction conditions
- **Better heat and mass transfer**
- **Improved safety** (smaller reactors, easier catalyst handling)
- Efficiency (smaller reactors; minimized downtime)

➤ Hydrogenation is one of the most popular reactions but so far it is mostly conducted in batch in the fine and specialty chemical industry




**Fig. 9.** Graphical comparison of the wall heat transfer intensity ( $UA_{th}/V_R$ ) with effective reactor volume ( $V_R$ ) for the 3 reactor technologies considered in this study. See Supporting information and Table S3 for details on the methodology

*Chemical Engineering Journal*  
441 (2022) 135951


# Oleic acid continuous hydrogenation

Batch : 1g Ni<sub>30%wt</sub>/Al<sub>2</sub>O<sub>3</sub> 120°C 20 bar 100 ml oleic acid  
Flow : 0.75g Ni<sub>30%wt</sub>/Al<sub>2</sub>O<sub>3</sub> 120°C 20 bar flow rate 0.05-1 ml.min<sup>-1</sup>  
Flow : 7.5g Pd<sub>2%wt</sub>/Al<sub>2</sub>O<sub>3</sub> (300 μm) 120°C 20 bar flow rate 10 ml.min<sup>-1</sup>



**K1 reactor single tube (12% capacity)**  
**Stearic acid (99.9%) production :**

**Conversion :** 99%  
**P** = 20 bars  
**ΔP** < 1 bar  
**ΔT** : isotherm (30 L.min<sup>-1</sup> circulating fluid)




**K5 capacity**  
**(Footprint : 0,4 x 0,4 x 0,4 meter)**

**Ni catalyst :** 1 kt/year  
**Pd Catalyst :** 6 kt/year

➤ This case study demonstrates the capability of KHIMOD equipment to intensify hydrogenation by a continuous process, delivering very high productivity in a small footprint

# Myrcene continuous hydrogenation to myrcane

## Myrcene



➤➤ Myrcene hydrogenation is a relevant case study to demonstrate how KHIMOD equipment can manage very exothermic hydrogenation reactions

A **challenging test reaction** for hydrogenation :

- **Very fast 1<sup>st</sup> step – slow 3<sup>rd</sup> step**
- **Highly exothermic**  
(  $\Delta_r H = 420 \text{ kJ.mol}^{-1}$  ;  $\Delta T_{\text{adiabatic}} = 1200\text{K}$  ; *nitrobenzene*  $\Delta_r H = 560 \text{ kJ.mol}^{-1}$   $\Delta T_{\text{adiabatic}} = 2500\text{K}$  )
- **With some selectivity issues** (dimerisation may occur to C<sub>20</sub>)
- Myrcene hydrogenation is similar to squalene hydrogenation, another case very relevant for the personal care industry


# Myrcene flow hydrogenation to alkanes : batch to flow

**Batch :** 0.5g Pd<sub>2%wt</sub>/Al<sub>2</sub>O<sub>3</sub> (50 μm)  
120°C 20 bar 100 ml myrcene 1M/heptane

**Flow :** 0.5g Pd<sub>2%wt</sub>/Al<sub>2</sub>O<sub>3</sub> (300 μm)  
120°C 20 bar, 8 ml/mn

Flow hydrogenation conducted with a very good control of the temperature

Conditions	Conc.	Yield	Kg/day/ g catalyst	Kg/day/ Lreactor	t/year on a K5
Batch	1 M	99,9 %	0,45	1 - 20	/
Flow	6 M (neat)	93 %	1,4	850	1 950
Flow	3 M	98 %	1,5	900	2 050



➤ This case study demonstrates the ability of KHIMOD reactors to manage very exothermic hydrogenation safely and open the possibility to reduce dramatically the use of solvents

# Case study : selective flow hydrogenation

Batch to flow challenge : Successive reactions :  $A \rightarrow B \rightarrow C$

Goal is to maximize B conversion and selectivity



►► Highly sensitive to flow issues (bypass – dead volume – diffusion)

➔ For this reaction flow cannot be as selective than model batch.

➔ But flow is the only industrial solution as the batch scale-up severely degrade performances

►► Selectivity vs Conversion curve compared with batch for flow validation of the reactor

## An example :



►► **Selective hydrogenation is one of the most challenging reaction as it require to design the equipment to deliver a very good plug flow behavior**

# Batch to flow hydrogenation in monotube reactor

## Reaction conditions :

- Catalyst Pd/C 5%wt - 200-300 $\mu$ m
- 0.5 M in cyclohexane
- 50°C – 5 bar
- Flow rate : 20 – 40 ml/min (catalyst content 3.5 g)

Contact time : Catalyst weight (g) / Liquid flow rate (ml.min<sup>-1</sup>)


- Batch : 1.25 g for 100ml,

Contact time : Catalyst weight (g) \* reaction time / Liquid volume

## Results :

- Very low pressure drop (below 1 bar)
- Very good pre-heating
- Very good agreement between ideal batch and flow for both conversion and selectivity

➤ This first experiment demonstrates that with a single channel set-up, the equipment has a very good plug flow behavior and can deliver both a high conversion and a high selectivity






# Lowering the fluid velocity impact the selectivity

## Reaction conditions

- Catalyst **Pd/C** 5%<sub>w</sub>t - 200-300 $\mu$ m
- 0.5 M in cyclohexane
- 50°C – 5 bar
- Flow rate : **3 – 40 ml/min** (catalyst content 3.5 g)

➔ 30 cm length reactive channel is perfectly suitable for fast reactions (contact time  $\tau < 30$ s)

➔ For slower reactions longer channels are required !



➔ This experiment shows that the very good plug flow behavior enabling to deliver a high selectivity and a high conversion is only reached when the flow rate is above 20 ml/mn

# Batch to continuous hydrogenation in a multitube reactor

## Reaction conditions


- Catalyst **Pd/Al<sub>2</sub>O<sub>3</sub>** 5%wt - 200-300µm
- # channels : 8 in parallel
- 0.5 M in cyclohexane
- 50°C – 5 bar
- Flow rate : 150 ml/min (ie. 19 ml/min per tube)
- catalyst content : 3.5 g per tube

**Pressure drop standard deviation : 2 %**

➔ **Very good agreement with ideal batch !**

➔ Very good translation from monotube as a function of flow rate

➔ Methodology usable for different catalysts with similar results



➔ This last experiment demonstrates how scale-up with the numbering-up approach can be achieved on selective hydrogenation thanks to the manifold developed by KHIMOD

# Continuous hydrogenation: cases studies

---

## KHIMOD

25 boulevard • Romain Rolland  
75014 Paris • France

